



Conjugation Protocol of Carboxylated Particles by Carbodiimide Chemistry

Protocol

The protocol is given for the coating of a fixed amount of 25 mg of particles with a special protein or antibody. It can be varied in the scale according to your individual requirements.

Material:

- 1. Particle suspension (surface: COOH or PEG-COOH) containing 25 mg of particles;
- 2. 0.5 M MES buffer (2-(4-morpholino) ethanesulphonic acid buffer), which was adjusted to pH 6.3 with 2.5 M Na2CO3:
- 3. 4 mg EDC (1-ethyl-3-(3-dimethylaminopropyl)-carbodiimide hydrochloride);
- 4. 8 mg NHS (N-hydroxysuccinimide);
- 5. 150 200 µg protein or antibody;
- 6. 0.01 M PBS buffer (pH = 7.4);
- 7. 200 µl of 25 mM glycine in 0.01 M PBS buffer.

Procedure:

- 1. Transfer the particle suspension into a 2 ml reaction tube;
- 2. Dissolve 4 mg EDC and 8 mg NHS in 0.5 M MES-buffer (pH = 6.3). The buffer volume should be a quarter of the initial volume of the particle suspension;
- 3. Incubate the suspension with continuous shaking for 45 min at room temperature;
- 4. Wash the activated particles by centrifugation (45,000 x g, 30 min) with 0.01 M PBS buffer (pH=7.4);
- 5. Add 0.01 M PBS buffer (pH = 7.4) containing 150 200 μ g protein or antibody;
- 6. Incubate the suspension with continuous mixing for 3 hours at room temperature;
- 7. Repeat step 4;
- 8. Add 200 µl of 25 mM glycine in 0.01 M PBS buffer incubate the suspension with continuous mixing for 30 min at room temperature;

CD Creative Diagnostics®



- 9. Repeat step 4;
- 10. Resuspend the particles in 1 ml PBS-buffer (pH = 7.4);
- 11. Stabilize the suspension by addition of 20 µl 1% sodium azide solution if necessary.

Note: This protocol is intended to provide general guidelines for the binding of biomolecules or related compounds. Further optimization may be required in order to achieve optimal functionality and stability from case to case.